

BodyFX™

by Inmode

BODY CONTOURING & CELLULITE TREATMENT




NON-INVASIVE SOLUTION TARGETING FAT WITHOUT LEAVING LAX SKIN. A.C.E. PROVIDES FULL HEAT CONTROL AND MAXIMUM SAFETY.



BodyFX treatments target problematic fatty tissues for ideal body contouring.

- **Radio-Frequency** provides precise and optimal heating of the skin for reduction, contraction and body shaping
- **Negative Pressure Massage** applies gentle vacuum pressure allowing for maximal depth treatment
- **Clinical studies** with histology's showing volume reduction and contraction



Controlled Heat

A.C.E. TECHNOLOGY BENEFITS

- Acquire. Control. Extend. Technology to safely optimize peak energy
- Pre-heating of tissue reduces threshold for selective apoptosis of fat cells
- Sophisticated tissue temperature and impedance monitoring achieves and maintains optimal temperature end points to prevent atrophic damage
- Active temperature monitoring and control with cut-offs for unprecedented safety with reduced risk of thermal injury


ROI Most other manufacturers sell body contouring as a stand alone \$80K+ system –BodyFX is one of your InMode modalities for an all-in-one aesthetic solution... with no consumables!

SPECIFICATIONS

Negative Pressure	Up to 500 mBar
RF Output Frequency	Up to 50 W
Output Frequency	1MHz
Temperature Cut-off	Real time temperature measurements set by operator up to 43°C




Dr. S. Mulholland - After 8 BodyFX Treatments



Dr. A. Carey - After 8 BodyFX Treatments



Dr. D. Duncan – 24 Year Old after 3 Months



Dr. S. Mulholland – Substantial Cellulite Improvement